



Recombinant Candida glabrata Cell wall synthesis protein KRE9 (KRE9)

Product Code	CSB-MP527945CZI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74683
Product Type	Recombinant Protein
Immunogen Species	Candida glabrata (strain ATCC 2001 / CBS 138 / JCM 3761 / NBRC 0622 / NRRL Y-65) (Yeast) (Torulopsis glabrata)
Purity	>85% (SDS-PAGE)
Sequence	DVNI VEPGPGDTFS GADGTVSINL KWVDNGAYPP LDKVAYYLFSLCYGPNTKIQ CVYKPDTKIT PDDLKDDAG TYSYKFDIQS SVVGSGQFYVQVFAWVDGQG YTLHYTPRIE LTNMGGPTTF ELLTYTDTIQ PVPQTSIRTG TDTNTQSFDA SSVYSLPYTE QTLKTRVAPI QQQPGTKVTA TTWSRRYPTS AVTFYSTFRS SLDQLSTLTP AWNYTVSSAV NYATPAAMPS DNGGWYNPIK RQSLSTRKLN VRYLSS
Source	Mammalian cell
Target Names	KRE9
Protein Names	Recommended name: Cell wall synthesis protein KRE9
Expression Region	17-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.