



Recombinant Schizosaccharomyces pombe GTPase-binding protein rid1 (rid1)

Product Code	CSB-MP527948SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74737
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNSLRSVKRL LGARSLSSDG LSESVNSSDR EDSLSFQTPS TPSEDEMPQC FEDLIKPQVL SSEASKSLLS LNEGKMEIV SSLQKQDQKK RNKLIWMRK KPCYEQPDDF ISYIVNTKID SIDESIIHKL ALLLRNEQVI WVEKFIHNGG FGVVFCTIDK ISRLEWREEL HDSILEQLLL CVKAMCTVQR GLECLSQSTY NVERLISLLL SKKQPCDFVV REVLVQIVHS FYKAHLDKKEE GAMQVFSLFS IKDDKDEEVE FLKNVRVDRP FRRWIELED VSRNVFWVWN HDSNVIDLEK QYNQVYEAPL GFIGGIETEA TSYVAAHIHL INELLEGLPL EKRQRKRELE QLSGLERIMG LRFRKSSQKF HKKLHEALRE WIVAAKIDDW PYKLVQTGST
Source	Mammalian cell
Target Names	rid1
Protein Names	Recommended name: GTPase-binding protein rid1
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.