



Recombinant Schizosaccharomyces pombe SCF E3 ubiquitin ligase complex F-box protein pof2 (pof2)

Product Code	CSB-EP527954SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74783
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MRVPNEVCFN ILSYLEADEL RCKSTVCTSW RNFIPTLWE KVVVFQNEAQL NNFFDTLQYS KDVSYYFRYL RKLNCSRVRK FLTDKHLMLM TLATGISRLN LSGCTRISEP LIGKLLYQNL NLVTINFSNI FSLPANILEY ISDNCPNLKA LNIGNCGLVE DTGMVQIIKR CPYLNRLIIP NCRKLTDVSL QILSEKEDLI ELDISGCEGF HNADTLRSLV SRNRGLKELS MDGCTELSHF ITFLNLNCEL DAMRALSLNN LPDLKDSIE LITCKFSKLN SLFLSKCIGL TDSSLLSLTK LSQSLTTLHL GHCEITDIG VQCLLKSKN ITYIDFGGCL RLSDIAVSAI AKLPYLQRVG LVKCICLTDL SVILLSGSFS RNLERVHLSY CIGLTAKSVS YLMYNCKTLK HLSVTGINSI LCTELRSFSR PIPDGINPSQ VPVFCAFTKV EIDLFREFIR NRI
Source	E.coli
Target Names	pof2
Protein Names	Recommended name: SCF E3 ubiquitin ligase complex F-box protein pof2 Alternative name(s): F-box and leucine-rich repeat protein pof2 F-box/LRR-repeat protein pof2 Short name= F-box protein pof2
Expression Region	1-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.