



# Recombinant Schizosaccharomyces pombe SCF E3 ubiquitin ligase complex F-box protein pof2 (pof2)

<b>Product Code</b>	CSB-YP527954SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O74783
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRVPNEVCFN ILSYLEADEL RCKSTVCTSW RNFIPTLWE KVVVFQNEAQL NNFFDTLQYS KDVSYYFRYL RKLNC SRVRK FLTDKHLMLM TLATGISRLN LSGCTRISEP LIGKLLYQNL NLVTINFSNI FSLPANILEY ISDNCPNLKA LNIGNCGLVE DTGMVQIIKR CPYLNRLIIP NCRKLTDVSL QILSEKEDLI ELDISGCEGF HNADTL SRLV SRNRGLKELS MDGCTELSHF ITFLNLNCEL DAMRALS LNN LPDLK DSDIE LITCKFSKLN SLFLSKCIGL TDSSLLSLTK LSQSLTTLHL GHCEITDIG VQCLLK SCKN ITYIDFGGCL RLS DIAVSAI AKLPYLQRVG LVKICLTDL SVILLSGSFS RNLERVHLSY CIGLTAKSVS YLMYNCKTLK HLSVTGINSI LCTELRSFSR PIPDGINPSQ VPVFCAFTKV EIDLFREFIR NRI
<b>Source</b>	Yeast
<b>Target Names</b>	pof2
<b>Protein Names</b>	Recommended name: SCF E3 ubiquitin ligase complex F-box protein pof2 Alternative name(s): F-box and leucine-rich repeat protein pof2 F-box/LRR-repeat protein pof2 Short name= F-box protein pof2
<b>Expression Region</b>	1-463
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.