



# Recombinant Schizosaccharomyces pombe Probable sphingolipid long chain base-responsive protein C736.15 (SPCC736.15)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP527976SXV  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | O74960   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MMNRAYSLRN SRAPTASQLL NPPPPSSSTR SGRLFAPLSH TMRINHAATF<br>TPDLAKRLAV LVKMEKNVMR SMEVTIRGRR DCARQLSYWG EDCDDDISDV<br>TDKLGVLFE MAELENYLID RYDQYRMTLK SIRNISSVQ PSREKKQKLL<br>DQIYALKHKD PESPRLVTME QELVREEAAC LVAEAQLSNT TREKFKQAMT<br>FNLDALHEHA EKLNLIATYG RHLLNLIDDT PVTPGEARPA YDGYETSRQI<br>VMDAEHALSS WVPSQPAVNF VKPQIEDDAQ SDARSWNEYE<br>AQGEPVPVEQ VTHLQDDVQS DTSEIENEP GMHPHVHAER MSRIASESSD<br>AVPQQTAVQV A |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | pil1   |
| <b>Protein Names</b>     | Recommended name: Probable sphingolipid long chain base-responsive protein C736.15<br>Alternative name(s): Protein kinase inhibitor C736.15  |
| <b>Expression Region</b> | 1-351  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.          |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |