



# Recombinant *Drosophila teissieri* Alpha-amylase-related protein (Amyrel)

<b>Product Code</b>	CSB-BP527989DMO
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O76260
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila teissieri</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	Q HNPWWGNRN TIVHLFEWKW SDIAQECESF LGPRGFAGVQ VSPVNIENIA AGRPWWERYQ PISYKLTTRS GNEEEFGDMV RRCNDVGVRI YVDVLLNHMS GDFDGVVVG T AGTEAEPREK SFPGVPTAQ DFHPTCEITD WNDRFQVQQC ELVGLKDLQ SSDWVRSKLI EFLDHLIELG VAGFRVDAAK HMASEDLEYI YSSLSNLNID HGFPHNSRPF IFQEVIDHGH ETVSRDEYKD LGAVTEFRFS EEIGNAFRGN NALKWLQSWG TGWGFLPSGQ ALTFVDNHDN QRDAGAVLNY KSPKQYKMAT AFHLAYPYGI SRVMSSFAFD DHDTPPPQDA QERIVSPEFD EDGACVNGWI CEHRWRQIYA MVGFKNAVRD TEVTGWWDNG DNQISFCRGN KGFLAINNNL YDLSQDLNTC LPEGTYCDVI SGSLIDGSCT GKSVTVNEYG FGYIHIGSDD FDGVLALHVD ARV
<b>Source</b>	Baculovirus
<b>Target Names</b>	Amyrel
<b>Protein Names</b>	Recommended name: Alpha-amylase-related protein EC= 3.2.1.1
<b>Expression Region</b>	20-493
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.