



Recombinant *Drosophila erecta* Alpha-amylase-related protein (Amyrel)

Product Code	CSB-YP527990DLK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O76265
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	Q HNPHWGGRN TIVHLFEWKW LDIAQECESF LAPRGFAGVQ VSPVNIIS AGRPWWERYQ PISYKLTTRS GNEEEFGEMV RRCNDVGVRI YVDVLLNHMS GDFDGVAVGT AGTQAEPRKK SFPGVPYSAQ DFHPTCEITD WDRFQVQQC ELVGLKDLNQ SSDWVRSKLI EFLDHLIELG VAGFRVDAAK HMASEDLEYI YSSLSNLNIA HGFPHNSRPF IFQEVIDHGH ETVSRDEYRD LGAVTEFRFS EEIGNAFRGN NALKWLQSWG TGWGFLPSGQ ALTFVDNHDN QRDAGAVLNY KSPKPYKMAT AFHLAYPYGI SRVMSSFAFD DHDTPPPQDA QERIISPEFD ADGACVNGWI CEHRWRQVYA MVGFKNAVRD TEITGWWDNG DSQISFCRGT KGFLALNNNL YDLSQDLNTC LPAGTYCDVI SGSLIDGSC GKSVTVNEHG YGYIHIGSDD FDGVLALHVD AKV
Source	Yeast
Target Names	Amyrel
Protein Names	Recommended name: Alpha-amylase-related protein EC= 3.2.1.1
Expression Region	20-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.