



Recombinant *Drosophila melanogaster* RILP-like protein homolog (CG11448)

Product Code	CSB-MP527996DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O76878
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MPGFHLNEMG EMVLDAIDDI GVVDVYDLAS DIGKEYERIM DRFGTDAVSG LMPKIINTLE LLEALATKNE RENATIQLER DKVAQLESEK LEKAEFRRRF DKELELIEEQ WRSETNELVD LVSSLQDENK RLVKQTQDLQ SSSAQSSGLG ASLTESIISM TNHELHSALS DTQVLQRLKE QIYKQRDELK HRERELQDKY SELEHLNIQA ERLKASERDT RRRHKLMQAQ VKTLCEERAD FLAQLQDQSR EINQLRKRLG LAEKENEDLV ASYDDGQNDP NRPRYTTREL KELISERDEL LTTIDTLNEQ LAELKPPSQA KGKRQRHFSS SDDSDEDDDDG HVADNDDDDDD EEEEAAEANE LEPPAAGETP PGHDAPVQGP LPYEPDDAPW KKSSSESGIRK FFRKLFS DPS DGSNTFPKRS LATLSKMALS ATPGSVSASA AAK
Source	Mammalian cell
Target Names	CG11448
Protein Names	Recommended name: RILP-like protein homolog
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.