



Recombinant Arabidopsis thaliana AP-4 complex subunit sigma (At2g19790)

Product Code	CSB-MP528083DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82201
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGIRFILMVN KQGQTRLAQY YEWLTLEERR ALEGEIVRKC LARNDQQCSF VEHRNYKIVY RRYASLFFMV GVDDDENELA ILEFIHLLVE TMDKHFGNVC ELDIMFHLEK AHFMLEEMVM NGCIVETSKA NILSPIQLMD KAH
Source	Mammalian cell
Target Names	At2g19790
Protein Names	Recommended name: AP-4 complex subunit sigma Alternative name(s): AP-4 adapter complex subunit sigma Adapter-related protein complex 4 subunit sigma Sigma subunit of AP-4 Sigma4-adaptin
Expression Region	1-143
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.