



# Recombinant Arabidopsis thaliana NPL4-like protein 2 (At2g47970)

<b>Product Code</b>	CSB-EP528087DOA-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O82264
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MMMLRIRSRD GLERVTAEGA HITVSQLKTL IADQLQIPLH KQTLSTNRDL LLAKTPADLL AFTDLTDPNL PLSSLNLGHG SMLYLAYDGE RSIPGAPPVT PAGSFGRKMT VDDLIARQMR VTRQETSHCD SVSFDRDAAN AFQHYVNESL AFAVKRGGFM YGTVTEEGQV EVDFIYEPQ QGTEANLILM RDADEEKRV AIAMGLGMRR VGFIFNQTVV QDKTEYTLN AEVLQAAELH AESELKEWVT AVVKLEVND GGADVFHFEAF QMSDMCIRLF KEEWFETEIM PDDDPKLSKM KKEVVVGVD LKEVDNDFFL VLVRILDHQQ PLSSTFPIEN RSSRATMRAL KTHLDRAKSL PLVKKMSDFH LLLFVAQFLD VSSDVPALAE CVRLQSPVPE GYALLIESMA NTC
<b>Source</b>	E.coli
<b>Target Names</b>	At2g47970
<b>Protein Names</b>	Recommended name: NPL4-like protein 2
<b>Expression Region</b>	1-413
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.