



Recombinant Salmonella typhimurium Probable rRNA maturation factor YbeY (ybeY)

Product Code	CSB-BP528337SXB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O87574
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	MSQVILDLQL ACENHAGLPD EAQFQRWLDG VIPQFQEEAE VTIRLVDEAE SHDLNLTYRG KDKPTNVLSF PFEAPAGIEM PLLGDLIICR QVVEQEAQEQ SKPLEAHWAH MVVHGSLLHLL GYDHIDDDEA EEMESLETEI MLAMGYEDPY IAEKIAE
Source	Baculovirus
Target Names	ybeY
Protein Names	Recommended name: Probable rRNA maturation factor YbeY
Expression Region	1-157
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.