



Recombinant *Schizosaccharomyces pombe* Uncharacterized transcriptional regulatory protein C2H10.01 (SPAC2H10.01)

Product Code	CSB-MP528431SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94392
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MTQKNTGPRR RTAVACDRRCR RRKIRCTGSD IPGQPCLACQ KAHADCHFST TSWRACRKKK APLHCPSTNG KTNVLPASYAS TPSLSPMTSS HFPYASDGTL ISATNHCDCC SYNPHYLPVN SNGSSPSHTS SLDSIPGYAS GSGPFKTPYR PVAVSESSSR SSFSMGSSSEF ASSTWSQSPM QSSLYVSPSS QPLDRFSPAS FPSKETLTSS SLSSSVPRSV SLSNFADSNT SNYPESSLLP AAKVGNNGFV GSNVLGMSY SNVEDSVRFQ TDSVSQMSSK SLEHLPMLSE MTLSSSASFG ASVSPKSTPG SNSTGAAVDV NSVHSNNGGSD LDYSSYLTDV SAVDSPSLDA DEVLRSFNLA NDPIGAVDPL TLEIDQGGSG QLLGQQSFAN DSLNAKQSSG KLGPLPPQSF NSSDYISFDE PSRYLSDDAS LWPDQYVDKV QPLATMKTCN PAVFSSNTSL DDMFFFIRDF DEDHPIQTAY
Source	Mammalian cell
Target Names	SPAC2H10.01
Protein Names	Recommended name: Uncharacterized transcriptional regulatory protein C2H10.01
Expression Region	1-480
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.