



# Recombinant Schizosaccharomyces pombe Probable phosphoglycerate mutase C1687.21 (SPAC1687.21, SPAC222.01)

<b>Product Code</b>	CSB-BP528449SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O94461
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKVFLIRHGQ TDQNKRGILQ GSVDTNLNET GRLQAKLLAQ RLLPLDIDQI FCSSMKRCRE TIAPYLELKP EVPIVYTDLI RERVYGDLEG MNVVEAKKLL NANHPDHYGE GLSHLTSRLL KFWDEYVPL QGKKKCVIVL CHGGVINVLR THFMEEKGFT FDKSKLESKI ITFNTSITEI EVTEHGTGHI HTFGDASHLE SEETGKLIA
<b>Source</b>	Baculovirus
<b>Target Names</b>	SPAC1687.21
<b>Protein Names</b>	Recommended name: Probable phosphoglycerate mutase C1687.21 EC=5.4.2.1 Alternative name(s): Phosphoglyceromutase C1687.21
<b>Expression Region</b>	1-209
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.