



Recombinant Schizosaccharomyces pombe 54S ribosomal protein L3, mitochondrial (mrpl3)

Product Code	CSB-YP528581SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O43042
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	SNSG KSLIENILSS KIAISKVYTL KKRIGLPESY PMKVLVQALI DKTVVQNPQY SNEKLWSIGK VFGDFYVLEY IKCKYPRLPE EAVHGMQNGW MGSKALSQLA TIWGLEVQRR EEGLQESFGK VLAKRADPEM LRKPGEIYKD EAARRAILAI LGGVYLHQGF NDTKKFINDH ILSKQLDPRT MLQLQWPRRQ LSRLCKRSL KEPVYRIIAE TGRKSREPVF VIGAYSGHHL LGQQQASSLN LAEQQAAYNS LISYYIHSPP FRQLPSDYLE VNDGKIPKAS EPYTSIPLVS HGQLVI
Source	Yeast
Target Names	mrpl3
Protein Names	Recommended name: 54S ribosomal protein L3, mitochondrial
Expression Region	27-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.