



Recombinant Schizosaccharomyces pombe Uncharacterized protein C947.09 (SPBC947.09)

Product Code	CSB-BP528585SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O43084
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPAKTRNVLV ACSDYYGPFY KDGENTGAFF LELLHPYLVF RDACFNVDIV TESGKIQFDD HSVAGPAIDK GSKGEEFLSY DDHIASGPEL SKAEKYVLEN KDDMFWRIVQ NSKTADEVNP DKYDIFFVAG GHATLFDLDFPK ATNLQKLGTS IYENGGVVAA VCHGPTLLPF MKRQTSDGSV SIVCGKDVTA FDRVAEDKSK LMEALKKYNL EVLDDMLNDA GANFIKSPNP FGDFVIADGR LVTGSNPASA TSTAKTALRV L
Source	Baculovirus
Target Names	hsp3103
Protein Names	Recommended name: Uncharacterized protein C947.09 EC= 3.2.-.-
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.