



Recombinant Putative T-box protein 33 (tbx-33)

Product Code	CSB-EP528600CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O45291
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MNYQNEQNAY QMQNQENVPP MSYGYPMP TP EQFQKMSYGQ NFQAPWDMQQ QIQMFQQYHQ NFQGGYSDPF TYQNQYAHQP KIVTVTLQDP QDLWKELHYL SNEQNVLNNG RKIFPALNYK VEHLNPESNY KVEILLRRMV PYQIQYSNGS WSRKNVQSKK TIAMKTEKVF VGEFTGQDIM RTGLDLSVK VFNIGSDNKK KVTPEEMSD AEKREYDTQY SKKKTSMLEV SNECKYIPVL IINEILPNQE LRLVGEFENE ITQFATVSSY KNHIVKALRT AANPTSRGDA KQEAQVQVKQ WRSSNKSLLE SLRPSMICTS VSTTAPSTNF HAPLQYPGTS SPSSNFAPMT PSTSFDSAYS SFNVTSSSTPE QMCYNPIPSM STDYSFCSD SATSSPPLQP TATSPEASQN QIKLEMNQYM
Source	E.coli
Target Names	tbx-33
Protein Names	Recommended name: Putative T-box protein 33
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.