



Recombinant Uncharacterized protein F14H3.2 (F14H3.2)

Product Code	CSB-EP528601CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O45363
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MTIPYLADLK DQTECKTCPK ILTRCKGSLY KVILHEFLMT AGAYFGVFLV FRFAINETQR EYAAEVFKKL KEQQNVCIPM QMMLAFFIAT VADQWEKIFE NVGYIENAAL AVATFLPDGK EKRDNNGNVI LAKVDNSNVR RNIIRYLVLS QILGIRDVSE LVKKRFANYD MIKATGVLQD HEEPLLKKVP CKTYAESFVP ITWIMSILQK FASKNEENLY YDTVYLEITD FYKKIIKLTR YDLIPIPLAY PQAVFLAVRI YFFFCLFTRQ HLDLEENWAL SHWGFPLLT LQFIFLVGCM KVAEILLNPM GQDDENFECN YVMDKNLFVG LTIVSSEHTE CPELEEVIGD DYVPWYPDDC KSKEEKNQEE LKKYLESVDF QAVTSSDQGE NDEVSTMMKV EQNSLLVCGR EKFYEGGGFR NSDVYPKQRA NYPH
Source	E.coli
Target Names	best-12
Protein Names	Recommended name: Uncharacterized protein F14H3.2
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.