



Recombinant Nuclear hormone receptor family member nhr-49 (nhr-49), partial

Product Code	CSB-EP528604CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O45666
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MDYFLDASKH IQLNQIDEES SDDFMDLVDP LAEPCAVCGD KSTGTHYGV I SCNGCKGFFR RTVLRDQKFT CRFNKRCVID KNFRCACRYC RFQKCVQVGM KREAIQFERD PVGSPTSGAS LNGTPFKKDR SPGYENGNSN GVGSNGMGQE NMRTVPQSSS VIDALMEMEA RVNQEMCNRY RRSQIFANGS GGSNGNDTDI QQGSDSGASA FAPPNRPCTT EVDLNEISRT TLLLMVEWAK TINPFMDLSM EDKIILLKNY APQHLILMPA FRSPDTTRVC LFNNTYMTRD NNTDLNGFAA FKTSNITPRV LDEIVWPMRQ LQMREQEFVC LKALAFHPE AKGLSNSSQI MIRDARNRVL KALYAFILDQ MPDDAPTRYG NILLAPALK ALTQLLIENM TLTKFFGLAE VDSLLSEFIL DDINDHSTAP VSLQQHLSSP TTLPTNGVSP LNPAGSVGSV SSVSGITPTG MLSATLAAPL AIHPLQSQDS ILNSEQNNHM L
Source	E.coli
Target Names	nhr-49
Protein Names	Recommended name: Nuclear hormone receptor family member nhr-49
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.