



# Recombinant Nuclear hormone receptor family member nhr-49 (nhr-49), partial

<b>Product Code</b>	CSB-MP528604CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O45666
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDYFLDASKH IQLNQIDEES SDDFMDLVDP LAEPCAVCGD KSTGTHYGV I SCNGCKGFFR RTVLRDQKFT CRFNKRCVID KNFRCACRYC RFQKCVQVGM KREAIQFERD PVGSPTSGAS LNGTPFKKDR SPGYENGNSN GVGSNGMGQE NMRTVPQSSS VIDALMEMEA RVNQEMCNRY RRSQIFANGS GGSNGNDTDI QQGSDSGASA FAPPNRPCTT EVDLNEISRT TLLLMVEWAK TINPFMDLSM EDKIILLKNY APQHLILMPA FRSPDTTRVC LFNNTYMTRD NNTDLNGFAA FKTSNITPRV LDEIVWPMRQ LQMREQEFVC LKALAFLHPE AKGLSNSSQI MIRDARNRVL KALYAFILDQ MPDDAPTRYG NILLAPALK ALTQLLIENM TLTKFFGLAE VDSLLSEFIL DDINDHSTAP VSLQQHLSSP TTLPTNGVSP LNPAGSVGSV SSVSGITPTG MLSATLAAPL AIHPLQSQDS ILNSEQNNHM L
<b>Source</b>	Mammalian cell
<b>Target Names</b>	nhr-49
<b>Protein Names</b>	Recommended name: Nuclear hormone receptor family member nhr-49
<b>Expression Region</b>	1-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.