



Recombinant Taraxacum officinale Root allergen protein

Product Code	CSB-BP528639TLH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O49065
Product Type	Recombinant Protein
Immunogen Species	Taraxacum officinale (Common dandelion) (Leontodon taraxacum)
Purity	>85% (SDS-PAGE)
Sequence	MAVAEFEITS SLSPSNIFKA FVIDFDTIAP KAEPETYKSI KTIEGDGGVG TIKSITYSDG VPFTSSKHKV DAIDSNNFSI SYTIFEGDVL MGIIESGTHH LKFLPSADGG SVYKHSMVFK CKGDAKLTDE NVSLMKEGLK KTFKAIETYV ISHPEAC
Source	Baculovirus
Protein Names	Recommended name: Root allergen protein Short name= RAP
Expression Region	1-157
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.