



Recombinant Mouse 5'-AMP-activated protein kinase subunit gamma-1 (Prkag1)

Product Code	CSB-MP528783MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O54950
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MESVAAESSP ALENEHFQET PESNNSVYTS FMKSHRCYDL IPTSSKLVVF DTSLQVKKAF FALVTNGVRA APLWDSKKQS FVGMLTITDF INILHRYYS ALVQIYELEE HKIETWREVV LQDSFKPLVC ISPNASLFDA VSSLIRNKIH RLPVIDPESG NTLYILTHKR ILKFLKLFIT EFPKPEFMSK SLQELQIGTY ANIAMVRTTT PVYVALGIFV QHRVSALPVV DEKGRVVDIY SKFDVINLAA EKTYNLDVS VTKALQHRSH YFEGVLKCYL HETLETIINR LVEAEVHRLV VVDEHDVVKG IVSLSDILQA LVLGTGGEKKP
Source	Mammalian cell
Target Names	Prkag1
Protein Names	Recommended name: 5'-AMP-activated protein kinase subunit gamma-1 Short name= AMPK gamma1 Short name= AMPK subunit gamma-1 Short name= AMPKg
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.