



Recombinant Arabidopsis thaliana Thioredoxin-like 1-1, chloroplastic (At1g08570)

Product Code	CSB-EP528948DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O64654
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	SSAITAQT TLRIGTAQKW WEKGLKDNMR EISSAQELVD SLTNAGDKLV VVDFDFSPGCG GCKALHPKIC QFAEMNPDVQ FLQVNYEEHK SMCYSLGVHV LPFFRFYRGS QGRVCSFSCT NATIKKFRDA LAKHGPDRCS LGPTKGLEEK ELVALAANKE LNFTYTPKPV PVEKEAATPD SNPSLPVPLP SMSSNDEKTL VSAGR
Source	E.coli
Target Names	At1g08570
Protein Names	Recommended name: Thioredoxin-like 1-1, chloroplastic Alternative name(s): Atypical cysteine/histidine-rich thioredoxin 4 Short name= AtAHT4 Lilium-type thioredoxin 1-1
Expression Region	73-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.