



Recombinant Schizosaccharomyces pombe

Chromatin modification-related protein png2 (png2)

Product Code	CSB-MP529198SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74736
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MGTSGIEIFA ALNDFTDAIV SVPESVCGKF TSLKEIDAQV RDIRQNVIQE IGVVLKNEKN DELSGEERCE RLQKTLKEIL PYSDSKICLA TDAMNNIKSC IDRLDAGFEY VELEIPQQLR LGYPDDRALM NYHSTVTPQT SERRRETRRH QNNQHSQQYS SQERSSSYNN FEDASSPQSS YHTPTKRRKN AVPRKSSSPP LSSTKHAPQS TERRPVRRSE SRLKQTNGEP LVKHDTLDSS DISREGEQLY CYCQQVSYGQ MIGCDNENCK REWFHLPCVG LVEPPKGIWY CKECEELAKS SESRQ
Source	Mammalian cell
Target Names	png2
Protein Names	Recommended name: Chromatin modification-related protein png2 Alternative name(s): ING1 homolog 2
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.