



Recombinant *Schizosaccharomyces pombe* Uncharacterized carboxypeptidase C337.07c (SPBC337.07c)

Product Code	CSB-YP529208SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74818
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	KS SRELGVKWSN DIRSWMTSYV GFNGETAAIS EEQLIWAEKI PDYNEEIVVR LHLEENNL DILQKAQSQN LDIWDYNFDH VDLRLKEENF DFWKSQYRSD ILINNL TETL FESIVPDTTN SPFSTEAFLQ AVENGHLNHE MFTSFTDIFF KSYQNLESIN SWLRLMASLY KDLSELPVG ITAEGRTLIG LKLNGRHPSD NGEKIRNKKV IIIQGGSHAR EWIGIPSVCY AAWQLLAKYD SDGHVRKLLD KFEWIFIPVL NVDGYEYTW S NDRLWSKNRQ PLNNSECFGI NLDANWAFGF NGNIDPCSNE YGGLSPFQAN ETMALFNLIT ESLSQEKKV VGFLDVHSYS QSVLWPYAYT CDLFPPDTEN FEELAIGLVK ELHRVNSRYY TYQQACIPYD GFHKHYLPGT AIDWVYFAAD VAWPFNIRLR DMGDYGYLLP AKQIVPTAKE FFAMILYGE FIAEYAF
Source	Yeast
Target Names	SPBC337.07c
Protein Names	Recommended name: Uncharacterized carboxypeptidase C337.07c EC=3.4.17.-
Expression Region	29-497
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.