



# Recombinant *Drosophila sechellia* Alpha-amylase-related protein (Amyrel)

<b>Product Code</b>	CSB-MP529233DMG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O76261
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila sechellia</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>Q HNPHWGGRN TIVHLFEWKW SDIAQECESF LGPRGFAGVQ  VSPVNIIS AGRPWWERYQ PISYKLTTRS GNEEEFGDMV RRCNDVGVRI  YVDVLLNHMS GDFDGVAVGT AGTEAEPRKK SFPGVPTAQ DFHPTCEITD  WNDRFQVQQC ELVGLKDLQ SSDWVRSKLI EFLDHLIELG VAGFRVDAAK  HMASEDLEYI YSSLSNLNID HGFPHNSRPF IFQEVIDHGH ETVSRDEYKD  LGAVTEFRFS EEIGNAFRGN NALKWLQSWG TGWGFLPSGQ ALTFVDNHDN  QRDAGAVLSY KSPRQYKMAT AFHLAYPYGI SRVMSSFAFD DHDTPPPQDA  QERIISPEFD EDGACVNGWI CEHRWRQIYA MVGFKNAVRD TEITGWWDNG  DNQIAFCRGN KGFLAINNNL YDLSQDLNTC LPAGTYCDVI SGSLIDGSC  GKSVTVNENG FGYIHIGSDD FDGVLALHVD AKV</p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Amyrel
<b>Protein Names</b>	Recommended name: Alpha-amylase-related protein EC= 3.2.1.1
<b>Expression Region</b>	20-493
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.