



Recombinant Cell death protein 6 (ced-6)

Product Code	CSB-EP529235CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O76337
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MAKDIYKTFK RSVSGIVGGN NINGEGSSSP STSAPQVKYR GGTGRTWIHP PDYLINGHVE YVARFLGCVE TPKANGSDVA REAIHAIRFQ RDLKRSEQTR ETAKLQKVEI RISIDNVIIA DIKTKAPMYT FPLGRISFCA DDKDDKRMFS FIARAEGASG KPSCYAFTSE KLAEDITLTI GEAFDLAYKR FLDKNRTSLE NQKQIYILKK KIVELETENQ VLIERLAEAL RANSKADYEN TGPPIYPGLG PPALPLSPMP QGPPPNIPPS SIYSMPRAND LPPTMAPTL PQISTSSNGA SPSVSPASTS PSGPAPSIPP PRPPALAPPP PVAPRRNPVV SPKNSTAGLL DGLELGS AEP AKKAPSNIFD DSFDPRAGEK KSTAAEYNPF GADFLSGIQN GKEAPPSASA ELLASEAIAR LPKPESSSVP PKKTAAEYDA MINEVEKKLA AMSSGSFEFG QLQTGDLGGI EGESDYGTPS DRLNPKMMNL KQ
Source	E.coli
Target Names	ced-6
Protein Names	Recommended name: Cell death protein 6 Alternative name(s): Candidate adapter protein ced-6
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.