



# Recombinant *Drosophila serrata* Alpha-amylase-related protein (Amyrel)

<b>Product Code</b>	CSB-MP529236DMH
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O76459
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila serrata</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QHNPQWWGNR NTIVHLFEWK WSDIAEECET FLAPRGFAGV QVSPVNENII SSGRPWWERY QPISYKLTR SGNEEFADM VRRCNVVGIR IYVDVLLNHM SGDFDGVAVG TAGTEAEPSK KSFPGVPYTA QDFHPSCEIT DWNDRFQVQE CELVGLKDLN QHSDYVRSKL IEFDLHLIEL GVAGFRVDAK KHMASEDLEY IYGSLSNLNI DHGFPHNARP FIFQEVIDHG HETVSREEYN ELGAVTEFRF SEEIGKAFRG NNALKWLQSW GTDWGFLKSE QALTFVDNHD NQRDQGAVLN YKSPRQYKMA TAFHLAYPYG ISRVMSFAF DDHDTPPPQD AQENIISPEF DEDGVCVNGW ICEHRWRQIY AMVGFKNAVR DTELHGWWDN GDNQISFCRG NKGFLAVNNN LYDLSQELNT CLPAGEYCDV TSGSLIDGAC TGKSVTVNEH GYGYIHIGSD DFDGVLALHV NAKV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Amyrel
<b>Protein Names</b>	Recommended name: Alpha-amylase-related protein EC= 3.2.1.1
<b>Expression Region</b>	21-494
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.