



Recombinant Enterobacteria phage RB69 DNA polymerase processivity component (45)

Product Code	CSB-EP529270ELM-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O80164
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage RB69 (Bacteriophage RB69)
Purity	>85% (SDS-PAGE)
Sequence	MKLSKDTIAI LKNFASINSG ILLSQGKFIM TRAVNGTTYA EANISDEIDF DVALYDLNSF LSILSLVSDD AEISMHTDGN IKIADTRSTV YWPAADKSTI VFPNKPIQFP VASVITEIKA EDLQQLLRVS RGLQIDTIAI TNKDGKIVIN GYNKVEDSGL TRPKYSLTLT DYDGSNNFNF VINMANMKIQ PGNKYKVMLWG AGDKVAAKFE SSQVSYVIAM EADSTHDF
Source	E.coli
Target Names	45
Protein Names	Recommended name: DNA polymerase processivity component Alternative name(s): DNA polymerase accessory protein 45 Gp45
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.