



Recombinant Arabidopsis thaliana Putative F-box protein At1g60370 (At1g60370)

Product Code	CSB-EP529288DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O80758
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MNGGEKLESI PIDLIEIHS RLPAESVARF RCVSKLWGSM FRRPYFTELF LTRSRARPRLL FVLQHNRRKW SFSVFSSPQN QNIYEKPSFV VADFHMKFSV STFPDFHSCS GLIHFSMMKG AYTVPVVCNP RTGQYAVLPK LTRTRYENSY SFGYDPIEK QIKVLFMSDP DSGDDHRILT LGTTEKMLGR KIECSLTHNI LSNEGVCING VLYYKASRIV ESSDDDTSD DDDDDHERSD VIVCFDFRCE KFEFIVICFY GQLINSVQLS LKQKSHQKLD LII
Source	E.coli
Target Names	At1g60370
Protein Names	Recommended name: Putative F-box protein At1g60370
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.