



Recombinant Solanum tuberosum Cysteine synthase, chloroplastic/chromoplastic

Product Code	CSB-EP529306FIG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O81155
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	AVSVPTKS STEIEGLNIA EDVTQLIGNT PMVYLNTIAK GCVANIAAKL EIMEPCCSVK DRIGFSMIVD AEEKGLISPG KTVLVEPTSG NTGIGLAFIA ASRGYKLILT MPASMSLERR VILKAFGAEL VLTDPKGMK GAVSKAEEIL NNTPDAYILQ QFDNPANPKI HYETTGPFIW EDTKGKIDIL VAGIGTGGTI TGTGRFLKEQ NPNIKIIGVE PTESNVLSSG KPGPHKIQGI GAGFIPGNLD QDVMDEVIEI SSDEAVETAR TLALQEGLLV GISSGAAALA AIQVGKRPEN AGKLGIVVFP SYGERYLSSI LFQSIREECE KMKPEL
Source	E.coli
Protein Names	Recommended name: Cysteine synthase, chloroplastic/chromoplastic EC= 2.5.1.47 Alternative name(s): CSase B Short name= CS-B O-acetylserine (thiol)-lyase O-acetylserine sulfhydrylase OAS-TL B
Expression Region	53-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.