



# Recombinant Arabidopsis thaliana Protein RKD5 (RKD5)

<b>Product Code</b>	CSB-MP529319DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O81791
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVDQGFFTLK KEKKKNILIK LFYVKILGFC FCDSWSSSDM AHSLTSLAVF QSVIRKEMVR SLHVESVEI EREFWFKSKS CYVEKKAKPL FRSEDFRRPE ISEGSVFGTW RCIFVFRFNH SLPRFPTLLC LSRNPKLEDI PNLANELKFI SELKPSKIYE EEQCSSSTEG YYNSDLPKPR KLVLKQDLNC LPDSETESEE SVNEKTEHSE FENDKTEQSE SDAKTEILKK KKRTPSRHVA ELSLEELSKY FDLTIVEASR NLKVGLTVLK KCREFGIPR WPHRRIKSLD CLIHDLQREA EKQQEKNEAA AMAVAKKQEK LETEKRNIK RPFMEIGIET KKFRQENFKK RHRASRAKKN QESLVTSSST
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RKD5
<b>Protein Names</b>	Recommended name: Protein RKD5 Short name= AtRKD5 Alternative name(s): RWP-RK domain-containing protein 5
<b>Expression Region</b>	1-370
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.