



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At2g29780 (At2g29780)

Product Code	CSB-EP529343DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82378
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAISETSDD GSHGGVPNKK PEELHKNPKE DDHQEEVEN HPPIPRQIPQ ALIRRTVALI KRCHYPSLSL LSKAFRIVIS SPELHQTRSS LNLTEPVLVA LIGFPPHSFP NWFILNHNIT RNIPLRLSAI GSLPPMNP GS AVVTIGYKMY VIGGLIGPNN PVKTVFVIDC RVHTCNYLPT MHRARYRAVA EVINGKIYVI GGCEKRYDDW IEVFDVVTGI WSTVPDRSHF MSSLPGGGFV TSVVMQNKIY ILDATCGLVY DPIDGTWESG ELGTTLKS YW YKPCCVIEDL LYSFDPYCLQ GSPINVDPN VMVWTPMMGT GIRAFDLDY FECKMANFGG KLMVFGATYN NPVTDTCWIE MALIKGETGP ILGMIDSMVP VLRSVYSPYI DLCRSVTF
Source	E.coli
Target Names	At2g29780
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At2g29780
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.