



Recombinant Arabidopsis thaliana Two-component response regulator ARR4 (ARR4)

Product Code	CSB-BP529355DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O82798
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MARDGGV SCL RRSEMMSVGG IGGIESAPLD LDEVHVLAVD DSLVDRIVIE RLLRIT SCKV TAVDSGWRAL EFLGLDNEKA SAEFDRLKVD LIITDYCMPG MTGYELLKKI KESSNFREVP VVIMSSENVL TRIDRCLEEG AQDFLLKPKV LADVKRLRSH LTKDVKLSNG NKRKLPEDSS SVNSSLPPPS PPLTISPESS PPLTVSTESS DSSPPLSPVE IFSTSP LSSP IDEEDDDVLT SSSEESPIRR QKMRSPGLD
Source	Baculovirus
Target Names	ARR4
Protein Names	Recommended name: Two-component response regulator ARR4 Alternative name(s): Response regulator 1
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.