



Recombinant Schizosaccharomyces pombe Meiotically up-regulated gene 167 protein (mug167)

Product Code	CSB-MP529668SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94427
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MVRYSLRKAH EYARDGENQE VLGEISQAIM YYGLAKSQFE KICQNTAQPL IKRCANNQIE ELMVRIRELR ESLPNKQTPI SMSMSTR LSP MYTSSLTPRF TSNSMVLPGD HNTTVRLKPS FREALMEDSE DDLYQMYSKF EMQVKKLSTN YGLAVAFSIV PMDDDQPSPS IFSSDESFLI VEGDDFSLGN SVEEDANATS KEDPAYQNTN EQIQPLSNFD ISQQEYLNNT RIPPYENEDLH LQHITQGTTD DNNVSKFLIP SYNDAKELSE EEMGRSHKRE ESFKRAFGHA SSSSESSIGEI TDSREDIQSN RLVNGSWENN DFTKEINNNF PDRSETPTLQ TIEAPTKLMK YTRRKS LFRF PFFRSISGKK KEEPME SGTD SLESSTAQIT VDSQLKIKQL ETQIATLQKQ LEQFQTSTLD QDLH
Source	Mammalian cell
Target Names	mug167
Protein Names	Recommended name: Meiotically up-regulated gene 167 protein
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.