



Recombinant *Schizosaccharomyces pombe* Chromatin structure-remodeling complex subunit rsc7 (rsc7)

Product Code	CSB-BP529689SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94522
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDSTAGSDAS TPTNHTVSSK KRRGRPPKGS YAAYGSSDDD DEEYSGRTRR RNRNTRPRVS APSSSSTVVP KDLYSHRATL EDDELNFGVV DPEGEKKVNE LGYLNGGREY RCRTFTCLGR GNRLYMLSTE PARAMGYRDS YLLFLKHRSL HKIIVDDSEK WDLIERNIIP HSYKGRAVGI VAARSIFREF GARIIVGGRR IVDDYWEGEF RARGFVEGEL ADPDDKLPPP GMPYNRNQYV AWHGASAVYH PQPSLEAQLP AAARKRKKEP PKDATWLFQH AKATAAYNND ITKYLQKQD IGYFEPHTNL LHVPLNTQPT KTHWIQAKTG TECPAPKLDI LVALNSDVAP QVSIANIPPS VYASCPLHVQ EAIRKRQEQE IRSIRLNSMY
Source	Baculovirus
Target Names	rsc7
Protein Names	Recommended name: Chromatin structure-remodeling complex subunit rsc7 Alternative name(s): Remodel the structure of chromatin complex subunit 7
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.