



Recombinant *Schizosaccharomyces pombe* tRNA-specific adenosine deaminase subunit tad2 (tad2)

Product Code	CSB-MP529715SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94642
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAGDSVKSAL IGIAGGPFSG KTQLCEQLLE RLKSSAPSTF SKLIHLTSFL YPNSVDRYAL SSYDIEAFKK VLSLISQGAE KICLPDGSCI KLPVDQNRIL LIEGYYLLLP ELLPYYTSKI FVYEDADTRL ERCVLQRVKA EKGDLTKVLN DFVTLSKPAY DSSIHPTREN ADIILPQKEN IDTALLFVSQ HLQDILAEMN KTSSSNTVKY DTQHETYMKL AHEILNLGPY FVIQPRSPGS CVFVYKGEVI GRGFNETNCS LSGIRHAELI AIEKILEHYP ASVFKETTLV VTVEPCLMCA AALKQLHIKA VYFGCGNDRF GGCGSVFSIN KDQSIDPSYP VYPGLFYSEA VMLMREFYVQ ENVKAPVPQS KKQRVLKREV KSLDLSRFK
Source	Mammalian cell
Target Names	tad2
Protein Names	Recommended name: tRNA-specific adenosine deaminase subunit tad2 EC=3.5.4.- Alternative name(s): tRNA-specific adenosine-34 deaminase subunit tad2
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.