



Recombinant Schizosaccharomyces pombe RNA polymerase II transcription factor B subunit 3 (pmh1)

Product Code	CSB-YP529724SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94684
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDDEGARKVD EKCPLCQADR YLNPNMKLLI NPECYHKMCE SCVDRIFFTG PAQCPTPGCN KILRKAKFRE QTFEDAQIER EVDVRKRISR IFNKGQQEFD SLQAYNDYLE EVEILTFNLI YKIDVEETEE KVKQYEQNR DSIAANSARA AAEARILAQN EILLKRQKQE AREAAIREHQ KEKERREQVE QQIIFDLATS GKDPNKIIQL SDSLKKQEN IASSVSNISR SSSILLSDVQ QVAEDTTPFS PLAGEKDGSK YFSYSKNTYQ DLYLEKVSHE PGRKCGFRIE DCHLRCLYEA FSGIDYDLES LKKLEVAS
Source	Yeast
Target Names	pmh1
Protein Names	Recommended name: RNA polymerase II transcription factor B subunit 3 Alternative name(s): CDK-activating kinase assembly factor MAT1 homolog RING finger protein pmh1 RNA polymerase II transcription factor B 38 kDa subunit RNA polymeras
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.