



Recombinant Schizosaccharomyces pombe Ubiquitin fusion degradation protein 1 (ufd1)

Product Code	CSB-MP529791SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42915
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MFGGSFFSSD DDGFSMMSQL RSAFHNNVNQ RFDTRYRCYP VAMIPGEERP NVNYGGKVL PPSALEKLSR LNVSYPLFD FENEAEEKK HGGVLEFIAE EGRVYLPYWM MTTLSLEPGD LVRVINTDIA QGSYVKLQPQ SVNFLDITDH RAVLENALRN FSTLTKSDIF EILYNDQVYQ IKVIDVQPDD SRHVSVVET DLVVDVDFPPI GYEESLQKNK QQNIAGVQGT MVTKIGYDEL VRQGDSNLMK GTGTKLNGKE VAEVPKINLL DVEKQCEPAP LILPLGTYFF GYPYKAPSIE EDSKKDPNLF KFEGAGTSLR ASRKTNGTMG KGSSDDPIDI DA
Source	Mammalian cell
Target Names	ufd1
Protein Names	Recommended name: Ubiquitin fusion degradation protein 1 Short name= UB fusion protein 1
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.