



Recombinant Schizosaccharomyces pombe 26S protease regulatory subunit 7 homolog (rpt1)

Product Code	CSB-YP529793SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42931
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPPKEDWEKY QKPVDTTEEN DKNPPPLDEG DIELLKSYAT GPYARELAAI KEDTEAVLKR INDTVGIKES DTGLAPISFW DVAADRQRMS EEQPLQVARC TKIINEQSA EKNAYVINLK QIAKFVVS LG ERVSPTDIEE GMRVGC DRNK YAIQLPLPPK IDPSVTMMQV EEKPDVTYGD VGGCKEQIER LREVELPLL SPERFVKLGI DPPKGIMLYG PPGTGKTLCA RAVANRTDAT FIRVIGSELV QKYVGEGARM VRELFEMART KKACIIFFDE IDAIGGARFD DGAGGDNEVQ RTMLELITQL DGFDPGRGNK VLFATNRPNT LDEALMRPGR IDRKVEFGLP DLEGRAHILR IHAKSMAIDK DIRWELIARL CPSQTGAELR SVCTEAGMFA IRARRRVATE KDFLDAVQKV VKGNQKFSST ADYMNMS
Source	Yeast
Target Names	rpt1
Protein Names	Recommended name: 26S protease regulatory subunit 7 homolog
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.