



Recombinant GDP-mannose 4,6 dehydratase 2 (gmd-2)

Product Code	CSB-YP529823CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O45583
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MEARNAEGLE SCIEKIQEVK LSSFAELKAF RERKVALITG ITGQDGSYLA ELLLSKGYKV HGIIRSSSF NTARIEHLYG NPVTHNGSAS FSLHYGDMTD SSCLIKLIST IEPTEIYHLA AQSHVKVSFD LPEYTAEVDA VGTLRLDDAI HACRLTEKVR FYQASTSELY GKVQEIPQSE LTPFYPRSPY AVAKMYGYWI VVNYREAYKM FACNGILFNH ESPRRGETFV TRKITRSVAK ISLRQQEHIE LGNLSALRDW GHAKEYVEAM WRILQQDTPD DFVIATGKQF SVREFCNLAF AEIGEQLVWE GEGVDEVGKN QDGVVRVKVS PKYYRPTEVE TLLGNPAKAR KTLGWEPKIT VPELVKEMVA SDIALMEADP MA
Source	Yeast
Target Names	gmd-2
Protein Names	Recommended name: GDP-mannose 4,6 dehydratase 2 EC= 4.2.1.47 Alternative name(s): GDP-D-mannose dehydratase Short name= GMD
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.