



Recombinant Drosophila melanogaster Accessory gland protein Acp62F (Acp62F)

Product Code	CSB-EP529828DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O46202
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MGWQGP KVDCTANGTQ TECPVACPET CEYSGNGPCV KMC GAPCVCK PGYVINERIP ACVLRSDCPK DVVRKEDMLL GVSNFKCF SR NYNCS
Source	E.coli
Target Names	Acp62F
Protein Names	Recommended name: Accessory gland protein Acp62F
Expression Region	25-115
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.