



Recombinant Arabidopsis thaliana Expansin-A4 (EXPA4)

Product Code	CSB-EP529856DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O48818
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	RIPGIYSGGA WQNAHATFYG GSDASGTMGG ACGYGNLYSQ GYGTNTAALS TALFNNGMSC GACFELKCAN DPQWCHSGSP SILITATNFC PPNLAQPSDN GGWCNPPREH FDLAMPVFLK IAQYRAGIVP VSYRRVPCRK RGGIRFTING HRYFNLVLIT NVAGAGDIVR ASVKGSRGTW MSLSRNWGQN WQSNVAVLVGQ ALSFRVTGSD RRTSTSWNMV PSNWQFGQTF VGKNFRV
Source	E.coli
Target Names	EXPA4
Protein Names	Recommended name: Expansin-A4 Short name= AtEXPA4 Alternative name(s): Alpha-expansin-4 Short name= At-EXP4 Short name= AtEx4 Ath-ExpAlpha-1.6
Expression Region	21-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.