



# Recombinant *Schizosaccharomyces pombe* Serine/threonine-protein kinase ppk25 (ppk25)

<b>Product Code</b>	CSB-BP530075SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O59697
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKPNTTNLRN ECWDTFSIPK RSQNIKINQS TKHQRSISDF VGTAGPGRQV GNWIIKKTIG AGSMGKVKLV VNILTGEKAA LKMIPFTPNN TSQTVRVQRE ALLGRLLRHP NICRVIDCIR TPACTYILFE YVPGGQLLEY ILARGKLD LARSFAMQLI NALVYLHKNF IVHRDLKIEN VLLTQDSRQV KLIDFGLSNF YSKDDLLRTY CGSLYFAAPE LLDKPYIGP EVDVWSLGVV IYVMVCGRVP FDDVSVPLMH SKIKSGKLEF PSYISEDCCS LIAAMLNVNP RKRCSELAQA KFPWLKNSF CLYLPIPLTS IPSTPSIRSH VFKPPFNKLV LQLLHEHGLA SIPELKHELY MAYIERKTTS LVCLYLLGVE SLAPALRIPT ALPPVYSRHQ RHHSEILGAM DLTEKITAMQ CPP
<b>Source</b>	Baculovirus
<b>Target Names</b>	ppk25
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase ppk25 EC= 2.7.11.1
<b>Expression Region</b>	1-423
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.