



Recombinant *Schizosaccharomyces pombe* Serine/threonine-protein kinase ppk25 (ppk25)

Product Code	CSB-EP530075SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O59697
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKPNTTNLRN ECWDTFSIPK RSQNIKINQS TKHQRSISDF VGTAGPGRQV GNWIIKKTIG AGSMGKVKLV VNILTGEKAA LKMIPFTPNN TSQTVRVQRE ALLGRLLRHP NICRVIDCIR TPACTYILFE YVPGGQLLEY ILARGKLD LARSFAMQLI NALVYLHKNF IVHRDLKIEN VLLTQDSRQV KLIDFGLSNF YSKDDLLRTY CGSLYFAAPE LLDKPYIGP EVDVWSLGVV IYVMVCGRVP FDDVSVPLMH SKIKSGKLEF PSYISEDCCS LIAAMLNVNP RKRCSELAQA KFPWLKNSF CLYLPIPLTS IPSTPSIRSH VFKPPFNKLV LQLLHEHGLA SIPELKHELY MAYIERKTTS LVCLYLLGVE SLAPALRIPT ALPPVYSRHQ RHHSEILGAM DLTEKITAMQ CPP
Source	<i>E. coli</i>
Target Names	ppk25
Protein Names	Recommended name: Serine/threonine-protein kinase ppk25 EC= 2.7.11.1
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.