



Recombinant *Schizosaccharomyces pombe* Histone deacetylase clr6 (clr6)

Product Code	CSB-BP530076SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O59702
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGFGKKKVSF FYDEDVGNVH YGPQHMKPH RVRMVHNLVV NYNLYEKLNV ITPVRATRND MTRCHTDEYI EFLWRVTPDT MEKFQPHQLK FNVGDDCPVF DGLYEFCSIS AGGSIGAAQE LNSGNAEIAI NWAGGLHHAK KREASGFCYV NDIALAALEL LKYHQRVLYI DIDVHHGDGV EEFFYTTDRV MTCSEHFKFGE YFPGTGHKID TGIGTGKNYA VNVPLRDGID DESYESVFKP VISHIMQWFR PEAVILQCGT DSLAGDRLGC FNLSMKGHSM CVDFVKSFNL PMICVGGGGY TVRNVARVWT YETGLLAGEE LDENLPYNDY LQYYGPDYKL NVLSNNMENH NTRQYLDSIT SEIENLRNL SFAPSVQMHK TPGDFTFENA EKQNIKKEEI MDERV
Source	Baculovirus
Target Names	clr6
Protein Names	Recommended name: Histone deacetylase clr6 EC= 3.5.1.98 Alternative name(s): Cryptic loci regulator 6
Expression Region	1-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.