



# Recombinant Schizosaccharomyces pombe Probable UTP--glucose-1-phosphate uridylyltransferase (SPCC794.10)

<b>Product Code</b>	CSB-BP530091SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O59819
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLHRRHFHFKS QSTLDFDSVA VSISASTMKN ELDKLVLSNR VSDKKTFGIQ MDNFFALYRR YLLHTVKGYE CDWDSIRPLG PEDMIDYGD LPLCKNAGKYL NRLAVVKLNG GMGNALGVNY PKAMIEVRDN QSFLDLSIRQ IEYLNRRYDV SVPFILMNSY DTNDETCVKL RKYAGCKIDI STFEQSRYP R VFVDSQLPVP KAAPSPIEEW YPPGHGDIFD ALVHSGTIER LLAQGKDYL F VSNIDNLGAS VDLNILSHVI DNQIEYSMEI TDKTKADIKV GILVNQDGLL RLLLETNQVPE QHREEFMSDK VFKYINTNNV WLYLPAVKRV VENRELNLDI MPNIETVYYN NEPARIIEFT TAIGSAISQF KKTEGIRVSR PRFISVKNSS DLFLVRCDL Y NVDHGSLKIE ESRLGFPPP V VRMSNEFKDI AELFCRIPYM PSMKDLVSL S ISGNVYFGRN VILKGNIVIV ASENTILCIP SNAVLENCVV TGNCKIMEC
<b>Source</b>	Baculovirus
<b>Target Names</b>	SPCC794.10
<b>Protein Names</b>	Recommended name: Probable UTP--glucose-1-phosphate uridylyltransferase EC= 2.7.7.9 Alternative name(s): UDP-glucose pyrophosphorylase Short name= UDPGP Short name= UGPase
<b>Expression Region</b>	1-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.