



Recombinant *Thermomyces lanuginosus* Lipase (LIP)

Product Code	CSB-MP530100TNC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O59952
Product Type	Recombinant Protein
Immunogen Species	<i>Thermomyces lanuginosus</i> (<i>Humicola lanuginosa</i>)
Purity	>85% (SDS-PAGE)
Sequence	EVSQDLFN QFNLFAQYSA AAYCGKNND A PAGTNITCTG NACPEVEKAD ATFLYSFEDS GVGDTVGF LA LDNTNKLIVL SFRGSRSIEN WIGNLNFDLK EINDICSGCR GHGFTSSWR SVADTLRQKV EDAVREHPDY RVVFTGHSLG GALATVAGAD LRGNGYDIDV FSYGAPRVGN RAFAEFLTVQ TGGTLYRITH TNDIVPRLPP REFGYSHSSP EYWIKSGTLV PVTRNDIVKI EGIDATGGNN QPNIPDIPAH LWYFGLIGTC L
Source	Mammalian cell
Target Names	LIP
Protein Names	Recommended name: Lipase EC= 3.1.1.3 Alternative name(s): Triacylglycerol lipase
Expression Region	23-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.