



Recombinant Putative ribosomal RNA methyltransferase F45G2.9 (F45G2.9)

Product Code	CSB-MP530133CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O62251
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MFSTKKSQGN LHKYIQRQST DEFAVKAREH NYRARSFAKL IEINEKFKFL KPESTVIDIG CAPGSWLQVV VQKCPNGYAS GVDLQNVLP I RGADILSLSD ITDPAVKLKI REKLAHRQVD VVLSDMAPNP TGDNATDHLR LIELCRSVFR LFSVENEIEL VKNGVYLCKI WDGSARAEFV RELSDRFSTV KTVKPTACRD NSAELYLFCR NFKK
Source	Mammalian cell
Target Names	F45G2.9
Protein Names	Recommended name: Putative ribosomal RNA methyltransferase F45G2.9 EC=2.1.1.- Alternative name(s): rRNA (uridine-2'-O-)-methyltransferase
Expression Region	1-214
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.