



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At4g39760 (At4g39760)

Product Code	CSB-EP530183DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O65663
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MISQEEVPQS TNHSLSFSSL PHEIVVSCLA RVSGSYYPKL CLVSKQFRSI ILSNEIYKAR SHLGTKENRL FVWLKLPTRS YPSWFALWIK PNETLTNDGP IKKQSTGNLL VPLPCSYNYQ VLVPSVIVGS ETYIVGGYDD ALSSSVWFYK NGKIHTLSKS PSMSVARIDA VVGQYPNIY VMGGCDSDES MNWGEVFNK TQTWEPLPDP GPEVRGQLVR KMKMQKKNVY VSSEKKDYIY DTEERTWKVT EAVFNFSWCV IEKVRYIYNN KNCWWLDTKS KDWRKIKGLD FLNKFRETDR IEIVNLDGKL VMIWDRFTLS KRNKKIWCAM IALEKCQGCE GIWGKIEWIG DVLMVPLSYS FLDCMVISI
Source	E.coli
Target Names	At4g39760
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At4g39760
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.