



Recombinant Arabidopsis thaliana Probable ribosome-binding factor A, chloroplastic (At4g34730)

Product Code	CSB-BP530184DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O65693
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	CMANPRRV KMKVAKQIMRE LSDMLLTDTV LQHAVLPEAA LGADRYLSSL TTISDVEVSN DLQVVKVYVS VFGDDRKDV AIAGLKSKAK YVRSELGKRM KLRLTPEVRF IEDEAMERGS RVIAILDKIK AEKSGGEGKT EPSDSTEDDQ DWEVDDPDED IIVVK
Source	Baculovirus
Target Names	At4g34730
Protein Names	Recommended name: Probable ribosome-binding factor A, chloroplastic
Expression Region	53-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.